



## Processing Item Information Sheet (PIIS)

### “Cell-Free DNA (cfDNA) Extraction from Whole Blood” Scheme [cfDNA23]

This sheet contains all the information on the **Whole Blood Processing Item** that you should be aware of to conduct the above-mentioned Scheme. **Please read carefully before performing any operation on the provided sample.**

### Processing Item Description

- Source material: Human blood from a healthy donor.
- Packaging: PAXgene ccfDNA tube.
- Date of Preparation: October 2023 – on the day of shipment.
- Testing of Biological Hazard: None. Blood was taken from a healthy donor, but it is recommended to consider the Item as potentially infectious.
- Biosafety level: All operations have been conducted in a BSL 2 environment.
- Homogeneity and Stability Information: Homogeneity of the Processing Item will be controlled in October 2023. Stability of the Processing Item is up to 10 days at room temperature (15–25°C) (data from PAXgene ccfDNA tube manufacturer).

### Instructions to Prepare the Processing Item for Extraction

- Any storage requirement between receipt and processing date: Store at **room temperature**. cfDNA extraction should be performed within 1 week of receipt.
- Requirements: Proceed to cfDNA extraction **immediately** after the whole blood centrifugation.

### Particular Handling/Safety Requirements

- Potential risks of Processing Item: All biological material should be manipulated as potentially hazardous.
- Individual protection equipment required: Standard laboratory equipment (laboratory coat, gloves).
- In case of puncture or cuts: Wash thoroughly with water and then disinfect during 10 minutes.
- In case of contact with the eye: Wash thoroughly with water or physiologic serum during 5 minutes.
- In case of contact with the mucus membranes and skin: Wash thoroughly with water.
- Measures to take in case of accidental spillage: Use disinfectant and thoroughly clean the effected surface.
- Waste elimination procedures: Waste generated by healthcare activities, to eliminate in incinerable plastic containers.

### Scheme Specifications

- Please extract **cfDNA** from the Processing Item following your **usual routine cfDNA extraction method**.
- You will be asked to report information under the following scheme: **Cell-Free DNA (cfDNA) Extraction from Whole Blood**.
- Please be ready to enter the following information:
  - Plasma isolation (centrifugation speed, duration and temperature, brake);
  - Kit used (QIAamp circulating Nucleic Acid kit (Diagnostic Sample Preparation, DSP, option) (QIAgen); QIAasymphony DSP Circulating DNA Kit/QIAamp DNA blood mini kit (QIAgen); QIAamp MinElute Kit (QIAgen); MagMAX™ Cell-Free DNA Isolation Kit (Applied Biosystems-ThermoFisher Scientific); MagNA Pure Compact



(MPC) Nucleic Acid Isolation Kit I (Roche Diagnostics); Maxwell® RSC (MR) ccfDNA Plasma Kit (Promega); NextPrep-Mag™ ccfDNA isolation kits (Perkin Elmer); NucleoSnap® DNA Plasma (Macherey-Nagel); FitAmp plasma/ serum DNA isolation kit (Epigentek); GenElute™ Plasma/Serum Cell-Free Circulating DNA Purification Midi Kit (Sigma); Quick-cfDNA Serum & Plasma Kit (Zymo), other);

- Extraction method (salting out, magnetic bead-based, silica membrane-based, unsure, other);
- Volume of plasma used for ccfDNA extraction (ml);
- Use of RNase;
- Elution buffer composition (elution buffer from kit, water, TE, AE, other) and elution volume (µl);
- Extraction equipment.

## What to Submit

- Once you have extracted DNA from the blood (according to your extraction method), you pipet all the extracted DNA in the provided labelled Matrix 0.5 ml. The tube is already labelled with your Laboratory Number and the name of the Scheme (cfDNA23). **Make sure to properly close the tube to avoid evaporation or leakage.**
- IBBL requires a **minimum of 10 µl** of extract to perform the downstream analyses planned.
- As soon as extracted, the DNA tube must be shipped to the following address, by using the courier of your choice:

IBBL, PT Programme  
Biorepository – Laura GEORGES/Margaux HENRY  
1B Rue Louis Rech  
L-3555 Dudelange  
LUXEMBOURG  
Phone: +352 26970-521  
Email: [biorepository@ibbl.lu](mailto:biorepository@ibbl.lu)

- The extracted DNA can be temporarily stored at -80°C before shipment. In that case, the extracted DNA must be shipped to IBBL on dry ice.
- Please note that IBBL cannot receive your sample on Saturdays nor Sundays.
- Your data must be submitted online on the PT website <http://biospecimenpt.ibbl.lu/> using the login information (Laboratory Number and Password) provided to you via email after the registration to the “Cell-Free DNA (cfDNA) Extraction from Whole Blood” Scheme.
- Please complete the questionnaire as accurately as possible, adding any relevant detail and comment in the appropriate section.

## Timelines

<i>Shipment of the extracted DNA to IBBL</i>	<i>Data Submission</i>	<i>Downstream analyses &amp; Report preparation</i>	<i>Reports available</i>
<b>Before</b> 10 NOV 2023	17 NOV 2023, <b>latest</b>	20 NOV 2023– 31 JAN 2024	29 FEB 2024

In case of doubts in the completion phase, please contact IBBL at [ISBERPT@ibbl.lu](mailto:ISBERPT@ibbl.lu)